

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: *H. pylori* Antigen, Unpurified

Catalog Number: HP-1001

Lot Number:

Quantity: _____ mg (= _____ ml)

Protein Concentration: 1.50 mg/ml by Pierce BCA protein assay

Product Form: Detergent extract of whole *H. pylori* cell sonicates. In 0.015M KPO₄ buffer, 0.85% NaCl, 0.1% NaN₃, pH = 7.40.

Specificity: *H. pylori* antibody reactive antigens

Storage Condition: -70°C or below

Cross-Reactivity: *H. pylori* specific. The antigens were extracted from a pure *H. pylori* culture derived from one single *H. pylori* colony.

Suggested Use: Microtiter Wells EIA or other immunoassays

Caution: The *H. pylori* antigen is an extract from whole *H. pylori* cells. This material should be handled as if it were capable of transmitting infection. Please use good laboratory practice techniques.