

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Anti- *H. pylori* Antibody, Affinity Purified

Catalog Number: HP-1005

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml by A_{280}

Form: *H. pylori* antigens-Affinity column chromatography purified immunoglobulin fractions in 0.015 M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN_3 .

Purity: > 95% by SDS-PAGE

Storage: At 2-8°C. DO NOT FREEZE.

Specificity: Reactive to *H. pylori* antigens.

Cross-Reactivity: *H. pylori* antigen specific. This antibody preparation is derived from goats immunized with antigens extracted from a pure *H. pylori* culture. It is immuno-reactive toward multiple *H. pylori* antigen as judged by western blot analysis and ELISA test.