

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Anti-Human IgM (μ)-HRP Conjugate

Catalog Number: MB 10123

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Form: Liquid

Buffer: 0.05M potassium phosphate buffer, 0.85% sodium chloride, with Preservative, pH=7.20

Storage: 2-8°C. DO NOT FREEZE

Suggested Use: Dilute this conjugate with BioCheck's Conjugate Diluent (P/N: MB 10149) to an appropriate dilution and use it as a secondary detection antibody.