

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CA 125 Standard/Sample Diluent

Catalog Number: MB-10025

Lot Number:

Quantity: _____ml

Product Form: In buffer solution containing 1.0% (v/v) ProClin-300 as preservative.

Storage: At 2-8°C. **DO NOT FREEZE.**

Expiration: 06/2011

Suggested Use: For use in the dilution of CA 125 Standards or patient samples in the CA 125 ELISA test.