

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** CA 15-3 Standard Diluent

**Catalog Number:** MB-10025

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In bovine serum containing 1.0% ProClin-300 as preservative.

**Storage:** At 2-8°C. **DO NOT FREEZE.**

**Expiration:**

**Suggested Use:** For use in the dilution of CA 15-3 Standards or patient samples in the CA 15-3 ELISA test.