

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Free β -hCG Standard/Sample Diluent

Catalog Number: MB-10026

Lot Number:

Quantity: _____ ml

Product Form: In buffer solution containing 1.0%(v/v) ProClin-300 as preservative.

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in the dilution of Free β -hCG Standards or patient samples in the Free β -hCG ELISA test.