

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Free β hCG Standard/Sample Diluent

Catalog Number: MB-10026

Lot Number:

Quantity: _____ ml

Product Form: In bovine serum containing 1.0% ProClin-300 as preservative.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of free beta-hCG Standards or patient samples in the hCG ELISA test.