

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** General Conjugate Diluent For  
Antibody-HRP Conjugate Dilution Use ( **$\beta$ -2 Microglobulin**)

**Catalog Number:** MB-10030

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.05M Tris buffer, pH=7.60, containing bovine serum proteins, preservatives and red dye.

**Storage:** At 2-8°C. **DO NOT FREEZE.**

**Suggested Use:** For use in the dilution of  $\beta$ -2 Microglobulin antibody-HRP conjugate in the  $\beta$ -2 Microglobulin ELISA assay.