

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Fatty Acid Binding Protein (FABP)
Antibody-HRP Conjugate Diluent

Catalog Number: MB-10030

Lot Number:

Quantity: _____ ml

Product Form: In 0.05M Tris buffer, pH=7.60, containing bovine serum proteins, preservatives and red dye. Filtered through a 0.2 μ membrane.

Storage: At 2-8°C. **DO NOT FREEZE.**

Expiration Date: _____.

Suggested Use: For use in the dilution of FABP antibody-HRP conjugate in the FABP ELISA assay.