

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Conjugate Diluent For Human Growth Hormone  
Antibody-HRP Conjugate Dilution Use (HGH)

**Catalog Number:** MB-10030

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.05M Tris buffer, pH=7.60, containing animal serum  
proteins, preservatives and red dye.

**Storage:** At 2-8°C. **DO NOT FREEZE.**

**Suggested Use:** For use in the dilution of HGH antibody-HRP conjugate in  
the HGH ELISA assay.