

# **BIOCHECK, INC.**

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## **PRODUCT SPECIFICATIONS**

**Product Name:** TMB Reagent, One Step

**Catalog Number:** MB-10051

**Lot Number:**

**Quantity:** \_\_\_\_\_ml

**Product Form:** Liquid containing H<sub>2</sub>O<sub>2</sub> and TMB.

**Storage:** At 2-8°C.

**Expiration:**

**Suggested Use:** For use as substrate/chromogen in HRP enzymatic reaction.