

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** T3 Antibody Diluent

**Catalog Number:** MB-10057

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.015 M KPO<sub>4</sub> Buffer, pH = 7.40, containing 0.25% (w/v) proteins, 0.5% (v/v) Proclin-300, and 0.02% (w/v) red dye.

**Storage:** At 2-8°C. DO NOT FREEZE.

**Suggested Use:** For use in the T3 antibody dilution in BioCheck's T3 EIA test.