

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: T-3/T-4 Conjugate Concentrate Buffer

Catalog Number: MB-10058

Lot Number:

Quantity: _____ ml

Product Form: In 0.05M Tris buffer, pH=7.60, containing bovine serum proteins and preservatives.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of T3-HRP or T4-HRP conjugate in the 11X concentrate form in the BioCheck's T3 or T4 EIA test.