

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Estradiol-HRP Conjugate Diluent

Catalog Number: MB-10058

Lot Number:

Quantity: _____ ml

Product Form: In 0.05M Tris buffer, pH=7.60, containing proteins and preservatives.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of Estradiol-HRP conjugate dilution in the BioCheck's Estradiol EIA test.