

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Estradiol-HRP Conjugate Diluent

**Catalog Number:** MB-10058

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.05M Tris buffer, pH=7.60, containing bovine serum proteins and preservatives.

**Storage:** At 2-8°C. **DO NOT FREEZE.**

**Suggested Use:** For use in the dilution of Estradiol-HRP conjugate dilution in the BioCheck's Estradiol EIA test.