

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Concentrated Wash Buffer (20x)

Catalog Number: MB-10073

Lot Number:

Quantity: _____ ml

Product Form: 100 mM Potassium Phosphate Buffer, 8.5% NaCl, 0.5% Tween 20, pH = 7.0

Storage: At 15-30°C. **DO NOT FREEZE.**

Suggested Use: Working Wash Buffer: Preparation of 1X Wash Buffer from 20X Stock. Add 50 ml of 20X Wash Buffer Stock to 950 ml of DI H₂O. The Working Wash Buffer is stable at 2-8°C for 30 days. NOTE: Any crystals that may be present due to high salt concentration must be redissolved at room temperature before making the dilution.