

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Conjugate Diluent For Antibody-HRP Conjugate Dilution Use

Catalog Number: MB-10076

Lot Number:

Quantity: _____ml

Product Form: Buffer solution, pH=6.20, containing proteins, preservatives and red dye.

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in the dilution of antibody-HRP conjugate in the ELISA assay.