

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Troponin I Standard Diluent

Catalog Number: MB-10077

Lot Number:

Quantity: _____ ml

Product Form: In human serum, containing 1.0% (v/v) ProClin as preservative.

Storage: At 2-8°C.

Suggested Use: For use in the dilution of Troponin I Standards or patient samples in the Troponin I ELISA test.