

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CK-MB Antibody-HRP Conjugate Diluent

Catalog Number: MB-10121

Lot Number:

Quantity: _____ ml

Product Form: Buffer, pH 6.2 containing bovine serum proteins, normal mouse IgG, normal goat IgG, preservatives, and red dye. Filtered through a 0.2 μ membrane.

Storage: At 2-8°C. DO NOT FREEZE.

Expiration Date: _____

Suggested Use: For use in the dilution of CK-MB antibody-HRP conjugate in the CK-MB ELISA assay.