

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: C-Reactive Protein (CRP) Antibody-HRP Conjugate Diluent

Catalog Number: MB-10143

Lot Number:

Quantity: _____ ml

Product Form: In 0.05 M Tris buffer, pH = 7.60, containing bovine serum proteins, normal mouse IgG, preservatives, and red dye. Filtered through a 0.2 μ membrane.

Storage: At 2-8°C. DO NOT FREEZE.

Expiration Date: _____

Suggested Use: For use in the dilution of CRP antibody-HRP conjugate in the CRP ELISA assay.