

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Conjugate Diluent for Antibody-HRP Conjugate Dilution Use

**Catalog Number:** MB-10149

**Lot Number:**

**Quantity:** \_\_\_\_\_ml

**Product Form:** Buffer solution, pH=6.70, containing proteins, preservatives and red dye.

**Storage:** At 2-8°C. **DO NOT FREEZE.**

**Suggested Use:** For use in the dilution of antibody-HRP conjugate in the ELISA assay.