

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Troponin I Antibody-HRP Conjugate Diluent, Solution B

Catalog Number: MB-10150

Lot Number:

Quantity: _____ ml

Product Form: In 0.05M KPO₄ buffer, pH=7.20, containing bovine serum proteins, 0.1% ProClin 300 and red dye.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of Troponin I antibody-HRP conjugate in the Troponin I ELISA assay.