

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Alkaline Phosphatase (ALP) Conjugate Diluent for Antibody-ALP Conjugate Dilution Use (1X)

**Catalog Number:** MB-10158

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.235M Borate buffer, pH=8.40, containing bovine serum Proteins and preservatives.

**Storage:** At 2-8°C. DO NOT FREEZE.

**Suggested Use:** For use in the dilution of antibody-ALP conjugate in the ELISA assay.