

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Alkaline Phosphatase (ALP) Conjugate Diluent for Antibody-ALP Conjugate Dilution Use (2X)

**Catalog Number:** MB-10159

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.470M Borate buffer, pH=8.40, containing bovine serum Proteins and preservatives.

**Storage:** At 18-30°C. Do not store at 2-8°C or freeze since precipitation may occur.

**Suggested Use:** For use in the dilution of antibody-ALP conjugate in the ELISA assay. **Dilute with equal volume of distilled water before use.**