

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** CK-MB ELISA Standard Set

**Catalog Number:** MB-30098

**Source:** Human Brain

**Lot Number:**

**CK-MB Concentration:** 0, 7.5, 15, 50, 100, and 200 ng/ml

**Quantity:** \_\_\_\_\_ set

**Product Form:** In bovine serum containing 1.0% (v/v) ProClin-300 as preservative.

**Storage:** -20°C or below.

**Viral Tests:** The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.