

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CK-MB ELISA Standard Set

Catalog Number: MB-30098

Source: Human Brain

Lot Number:

CK-MB Concentration: 0, 7.5, 15, 50, 100, and 200 ng/ml, liquid frozen.

Quantity: _____ sets

Product Form: In buffer solution containing 1.0% (v/v) ProClin-300 as preservative.

Storage: -20°C or below.

Viral Tests: The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.