

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** CA 15-3 MoAb-HRP Conjugate Concentrate (22X)

**Catalog Number:** MB-40008C

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.05 M KPO<sub>4</sub> buffer, 0.85% NaCl, pH=7.20,  
containing protein stabilizers and preservatives.

**Storage:** At 2-8°C. Do not freeze.

**Expiration:**

**Dilution Factor:** 1:22 in the **CA 15-3 ELISA assay**