

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** T-3-HRP Conjugate Diluent

**Catalog Number:** MB-50003

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In 0.10M Tris buffer, pH=7.60, containing bovine serum proteins, ANS, preservatives and blue dye.

**Storage:** At 2-8°C. DO NOT FREEZE.

**Suggested Use:** For use in the dilution of 11 X T3-HRP conjugate in the T3 EIA assay.