

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CA 19-9 Assay Buffer

Catalog Number: MB-50009

Lot Number:

Quantity: _____ ml

Product Form: 0.015 M KPO₄ Buffer, 0.85% (w/v) NaCl, 0.01% (w/v) Mouse IgG, 0.1% (w/v) NaN₃, 0.25% (w/v) BSA, 0.01% (w/v) EDTA(Na)₂, 0.005%(w/v) green dye, pH = 7.40

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in CA 19-9 ELISA as Assay Buffer.