

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Progesterone-HRP Conjugate Diluent

Catalog Number: MB-50048

Lot Number:

Quantity: _____ ml

Product Form: In 0.10M Tris buffer, pH=7.60, containing bovine serum proteins, ANS, preservatives and 0.005% (w/v) blue dye.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of Progesterone-HRP Conjugate in the BioCheck's progesterone EIA assay.