

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CRP Standard/Sample Diluent

Catalog Number: MB-50059

Lot Number:

Quantity: _____ml

Product Form: In 0.015 M KPO₄, 0.85% NaCl, 0.5% (v/v) Proclin 300,
0.0025 % (w/v) yellow dye, 0.0025% (w/v) red dye, pH = 7.40

Storage: At 2-8°C