

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CRP Standard / Sample Diluent

Catalog Number: MB-50059

Lot Number:

Quantity: _____ ml

Product Form: [6% BSA in 0.015 M KPO₄, 0.85% NaCl, 0.5% (v/v) Proclin 300, 0.0025% (w/v) Yellow Dye, 0.0025% (w/v), Red Dye, pH = 7.40]

Storage: At 2-8°C.