

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Zero Buffer for AFP/IgE ELISA

**Catalog Number:** MB-60005/MB-60008

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** 0.05 M Tris Buffer, 2.0% (w/v) proteins, 1% (v/v) Proclin-300, 0.005% (w/v) yellow dye, pH = 7.60

**Storage:** At 2-8°C. DO NOT FREEZE.

**Suggested Use:** For use in the AFP/IgE ELISA as Zero Buffer