

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Zero Buffer for AFP/IgE ELISA

Catalog Number: MB-60005/MB-60008

Lot Number:

Quantity: _____ ml

Product Form: 0.05 M Tris Buffer, 2.0% BSA, 1% Proclin-300, 0.005% Yellow Dye,
pH = 7.60

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the AFP/IgE ELISA as Zero Buffer