

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Zero Buffer for IgE ELISA

Catalog Number: MB-60008

Lot Number:

Quantity: _____ ml

Expiration: 05/2011

Product Form: 0.05 M Tris Buffer, 2.0% (w/v) BSA, 1% (v/v) Proclin-300, 0.005% (w/v) Yellow Dye, pH = 7.60

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in the IgE ELISA as Zero Buffer.