

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: PSA Zero Buffer

Catalog Number: MB-60010

Lot Number:

Quantity: _____ ml

Product Form: In bovine serum containing 1.0% ProClin-300 as preservative, and 0.005% Blue Dye.

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in PSA, Free PSA, hCG, and Free β -hCG ELISA as Zero Buffer.