

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Free PSA Zero Buffer

**Catalog Number:** MB-60011

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In buffer solution containing 1.0% (v/v) ProClin-300 as preservative, and 0.005% (w/v) blue dye.

**Storage:** At 2-8°C. DO NOT FREEZE.

**Suggested Use:** For use in PSA, Free PSA, hCG, and Free  $\beta$ -hCG ELISA as Zero Buffer.