

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Free β -hCG Zero Buffer

Catalog Number: MB-60014

Lot Number:

Quantity: _____ ml

Product Form: In buffer solution containing 1.0% (v/v) ProClin-300 as preservative, and 0.005% (w/v) blue dye.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in PSA, Free PSA, hCG, and Free β -hCG ELISA as Zero Buffer.