

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: β -2 Microglobulin Standard/Sample Diluent

Catalog Number: MB-60027

Lot Number:

Quantity: _____ ml

Product Form: In 0.05 M Tris buffer, 3.0% BSA (w/v), 0.1% NaN₃ (w/v), and 0.005% (w/v) Yellow Dye, pH = 7.60.

Storage: At 2-8°C. DO NOT FREEZE.

Suggested Use: For use in the dilution of β -2 Microglobulin Standards or patient samples in the β -2 Microglobulin ELISA test.