

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Streptavidin-HRP Conjugate Diluent

Catalog Number: MB-60116

Lot Number:

Quantity: _____ ml

Storage: At 2-8°C. DO NOT FREEZE.

Expiration: _____.

Suggested Use: For use in the dilution of Streptavidin-HRP Conjugate in the enzyme immunoassay.