

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: 650-573-1968 Fax: 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Rabbit Monoclonal Anti-Mouse/Human Id4. Produced by Recombinant cloning of Antibody Gene (cDNA), Clone #82-12
Catalog Number:	BCH-9/#82-12-cDNA
Antigen Source:	Full Length Human Id4-GST
Lot Number:	
Quantity:	
Antibody Concentration (A_{280nm}):	0.05 mg/ml
Product Form:	DEAE Column Purified rabbit monoclonal anti- mouse/human Id4, in 0.015M KPO_4 buffer, 0.85% (w/v) NaCl, 0.1% (w/v) NaN_3 0.25% (w/v) BSA, pH 7.4.
Specificity:	Mouse Id4 and Human Id4
Cross Reactivity:	No cross-reactivity with Mouse/Human Id1, Mouse/Human Id2, Mouse/Human Id3, or Carrier Protein GST
Storage Condition:	At 2-8° C. <u>DO NOT FREEZE.</u>
Suggested Use:	For Western Blot and Immunohistochemistry use.
Suggested Dilution:	For Western Blot: 1:500. For IHC staining: 1:100. Optimal conditions for Western Blot and IHC should be determined by the individual researcher.