

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Cytomegalovirus IgM Antigen
Catalog Number:	70168
Lot Number:	
Protein Concentration:	
Quantity:	
Strain:	AD-169
Produced in:	HF Cell line
Product Form:	Liquid. In 0.1 M Glycine Buffer, pH = 9.3 to 9.7
Preparation:	This antigen is prepared by the purification of virus from human fibroblast cells. Purity is greater than 90%.
Infectivity:	Negative.
Storage Condition:	-70° or below
Applications:	For research use only.
Animal Source Products / Country of Origin:	Product manufactured in the United States using only domestic animal materials processed in U.S. facilities.
Hazard/Biohazard:	The antigen must be considered potentially infectious and handled accordingly.