

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: IgE MoAb HRP Conjugate

Catalog Number: 70188-HRP

Lot Number:

Quantity: _____ml

Antibody Concentration (A_{280nm}): _____mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH = 7.20, containing proteins, and preservatives.

Storage: 2-8°C. DO NOT FREEZE

Dilution Factor: Greater than 1:1,000 in the human IgE ELISA assay.