

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Human Chorionic Gonadotropin Antigen,
Highly Purified

Catalog Number: 70264 (4002)

Lot Number:

Source: Human Pregnancy Urine.

Quantity: _____ mg (= _____ ml)

Specific Immunological Activity: _____ IU/mg/ml

Reference Preparation: WHO, 3rd IS (75/537)

Purity: > 98%. As determined by SDS-PAGE.

Product Form: In 0.015M KPO₄ buffer, 0.85% NaCl,
0.1% NaN₃, pH=7.40

Level of Specific Contaminants: hCG 100%
hLH < 0.01%
hFSH < 0.01%
hTSH < 0.01%

Storage: - 20°C or below.