

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	<b>Rabbit Polyclonal Anti-Human Alpha Fetoprotein (AFP) (Total IgG Fractions)</b>
<b>Catalog Number:</b>	70269 (# 08)
<b>Antigen Source:</b>	Human Cord Blood AFP
<b>Lot Number:</b>	
<b>Quantity:</b>	_____ mg (= _____ ml)
<b>Antibody Concentration (A<sub>280nm</sub>):</b>	_____ mg/ml
<b>Product Form:</b>	DEAE column chromatography purified immunoglobulin fractions in 0.015M KPO <sub>4</sub> buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN <sub>3</sub> .
<b>Specificity:</b>	Human Alpha Fetoprotein (AFP)
<b>Storage:</b>	At 2-8° C. <u>DO NOT FREEZE</u> . Precipitation may occur upon freezing.
<b>Cross-Reactivity:</b>	No cross-reactivity with Human Albumin, CEA, Ferritin, hCG, PAP, and PSA.
<b>Suggested Use:</b>	For solid phase immobilization (capture) in the human AFP sandwich immunoassay.
<b>Pairing Information:</b>	This antibody pairs with AFP MoAb (P/N: 70448) in the human AFP sandwich immunoassay.