

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Sheep Polyclonal Anti-Human β-2 Microglobulin, Affinity Purified
Catalog No:	70320
Lot Number:	
Immunogen:	Highly purified human β -2 Microglobulin
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Product Form:	Affinity column chromatography purified sheep polyclonal anti-human β -2 microglobulin in 10 mM sodium phosphate buffer, pH = 7.20, containing 0.85% NaCl and 0.1% NaN ₃ .
Storage:	2-8° C
Specificity:	This antibody reacts with human β -2 Microglobulin.
Suggested Use:	For Enzyme Conjugation (Detection)
Pairing Information:	This antibody pairs with β -2 microglobulin MoAb (P/N: 70191) in the human β -2 microglobulin sandwich immunoassay.