

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human Beta-FSH
Catalog Number:	70325 (G-3001)
Lot Number:	
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Subclass:	IgG ₁
Product Form:	DEAE column chromatography purified immunoglobulin fractions in 0.015 M KPO ₄ buffer, pH = 7.40 containing 0.85% NaCl, 0.1% NaN ₃ .
Specificity:	Human β eta-FSH subunit
Affinity Constant:	1×10^{10} l/mol.
Storage:	2-8° C. <u>DO NOT FREEZE</u> . Precipitation may occur upon freezing.
Cross-Reactivity:	FSH 100%, Beta hCG < 1%, ihCG < 1%, LH < 1%, TSH < 1%
Suggested Use:	For enzyme conjugation (detection).
Pairing Information:	This antibody pairs with Alpha-FSH MoAb (P/N: 70327) in the human FSH sandwich immunoassay.