

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Polyclonal Anti Human-Alpha-TSH, Affinity Purified

**Catalog Number:** 70343 (G-7005)

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg ( = \_\_\_\_\_ ml )

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** TSH and hCG polyclonal affinity columns purified Immunoglobulins from goat anti-human TSH serum in 0.015 M  $KPO_4$  buffer, pH=7.40, containing 0.85% NaCl, 0.1%  $NaN_3$ .

**Purity:** Greater than 98% by protein gel electrophoresis

**Storage:** At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

**Specificity:** Human Alpha-TSH

**Cross-Reactivity:** Cross react with Alpha hCG, Alpha-LH, and Alpha-FSH.

**Suggested Use:** For enzyme conjugation (detection) or solid phase coating (capture).