

BIOCHECK, INC.
425 Eccles Ave
South San Francisco, CA 94080
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Mouse IgG (H+L),
Affinity Purified

Catalog Number: 70346 (3001)

Lot Number:

Quantity: _____ mgs (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Form: Mouse IgG affinity column purified goat polyclonal anti-
mouse IgG(H+L). In 0.10 M potassium phosphate buffer,
pH = 7.80 containing 0.05% NaN₃.

Purity: Greater than 98% pure

Storage: At 2-8°C. DO NOT FREEZE.

Specificity: This antibody reacts with heavy chains on mouse IgG,
and with light chains on all mouse immunoglobulins. No
reactivity is observed to non-immunoglobulin mouse
serum proteins.

Note: Slight precipitation of antibody may occur on storage at
2-8°C. Refilter if necessary.