

# BIOCHECK, INC.

425 Eccles Ave  
South San Francisco, CA 94080  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Polyclonal Anti-Rabbit IgG (H+L),  
Affinity Purified

**Catalog Number:** 70347

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Form:** Rabbit IgG affinity column purified goat polyclonal anti-rabbit IgG (H+L). In 0.10 M potassium phosphate buffer, pH = 7.80 containing 0.05% NaN<sub>3</sub>.

**Purity:** Greater than 98% pure

**Storage:** 2-8°C, DO NOT FREEZE

**Specificity:** This antibody reacts with heavy chains on rabbit IgG, and with light chains on all rabbit immunoglobulins. No reactivity is observed to non-immunoglobulin rabbit serum proteins.

**Note:** Slight precipitation of antibody may occur on storage at 2-8°C. Centrifuge or refilter if necessary.